

Marker Kit

Gold fiducial markers


Accurate localization for soft tissue targets



Sterile Kit with 5 Markers (Loose)



Sterile Kit of 4 Needles with Marker in Each



Sterile Kit with Single Marker in Needle



Best Medical's gold markers are available in a variety of diameters and lengths to meet your needs for soft tissue localization such as prostate, liver, pancreas, lung, and breast.

- Best® gold markers are available now in tissue-absorbable strands (patent pending)
- Standard strand consists of two markers at 2 cm apart center to center
- Special marker strands with custom spacing are also available

Diameters	Varying Lengths
0.8 mm	1.5 mm
	3 mm
	5 mm
1.0 mm	7 mm
	1.5 mm
	3 mm
1.2 mm	5 mm
	7 mm
	3 mm
	5 mm

Best® Fiducial Marker Kit

Best® Gold Markers




- Excellent visualization under X-ray, Fluoro, CT, MRI, Ultrasound and other imaging modalities
- Sold individually, stranded, and in pre-loaded, pre-sterilized needle sets
- Can be used for brachytherapy, IMRT and IGRT

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
 phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



Anderson Marker Kit

Part # 1074-0002



Kit contains a 25 cm long implanter (designed to limit implantation to 1cm) and 10 stainless steel markers. The markers are easily identified on visualization films and also aid in treatment planning.

Markers—Part # 179-1: additional markers sold by vials (each vial contains 50 markers)

Adjustable Marker Kit

Part # 1074-0011



Kit contains a 20 cm long implanter with adjustable stop, one spare needle, one Allen wrench and 10 stainless steel markers.

Markers—Part # 179-1: additional markers sold by vials (each vial contains 50 markers)



NOTE: Gold markers (Part # 351) are also available to use with Anderson Marker Kit or Adjustable Marker Kit

CAUTION: This product should be carefully inspected by your staff prior to clinical use and should be used by authorized personnel only.

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com


Best[®]
healthcare for everyone

Sterile Gold Marker Kits and Radiopaque Strands

Part Number	Description
Sterile Gold Markers With Forceps	
782-1	0.8 mm diameter, 3 mm long (5 markers)
782-2	0.8 mm diameter, 5 mm long (5 markers)
782-3	0.8 mm diameter, 7 mm long (5 markers)
782-4	0.8 mm diameter, 10 mm long (5 markers)
784-1	1.0 mm diameter, 3 mm long (5 markers)
784-2	1.0 mm diameter, 5 mm long (5 markers)
784-3	1.0 mm diameter, 7 mm long (5 markers)
784-4	1.0 mm diameter, 10 mm long (5 markers)
785-1	1.2 mm diameter, 3 mm long (5 markers)
785-2	1.2 mm diameter, 5 mm long (5 markers)
785-3	1.2 mm diameter, 7 mm long (5 markers)
785-4	1.2 mm diameter, 10 mm long (5 markers)

Part Number	Description
Sterile Fiducial Marker Kit (1 per pack), 18 gauge, 10 cm Long Needle (0.8 mm OD)	
1323-1	18 gauge, 10 cm long Pre-waxed needle (1) preloaded with 0.8 mm OD, 3 mm long gold marker (1)
1323-2	18 gauge, 10 cm long Pre-waxed needle (1) preloaded with 0.8 mm OD, 5 mm long gold marker (1)
1323-3	18 gauge, 10 cm long Pre-waxed needle (1) preloaded with 0.8 mm OD, 7 mm long gold marker (1)
1323-4	18 gauge, 10 cm long Pre-waxed needle (1) preloaded with 0.8 mm OD, 10 mm long gold marker (1)
1323-5	18 gauge, 10 cm long Pre-waxed needle (1) preloaded with 0.8 mm OD, 1.5 mm long gold marker (1)
Sterile Fiducial Marker Kit (2 per pack), 18 gauge, 10cm Long Needle (0.8 mm OD)	
1428-1	18 gauge, 10 cm long Pre-waxed needle (2) preloaded with 0.8 mm OD, 3 mm long gold marker (1)
1428-2	18 gauge, 10 cm long Pre-waxed needle (2) preloaded with 0.8 mm OD, 5 mm long gold marker (1)
1428-3	18 gauge, 10 cm long Pre-waxed needle (2) preloaded with 0.8 mm OD, 7 mm long gold marker (1)
1428-4	18 gauge, 10 cm long Pre-waxed needle (2) preloaded with 0.8 mm OD, 10 mm long gold marker (1)
1428-5	18 gauge, 10 cm long Pre-waxed needle (2) preloaded with 0.8 mm OD, 1.5 mm long gold marker (1)

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
 phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



Sterile Gold Marker Kits and Radiopaque Strands—continued

Part Number	Description
	Sterile Fiducial Marker Kit (3 per pack), 18 gauge, 10 cm Long Needle (0.8 mm OD)
1325-1	18 gauge, 10 cm long Pre-waxed needles (3) preloaded with 0.8mm OD, 3 mm long gold marker (1)
1325-2	18 gauge, 10 cm long Pre-waxed needles (3) preloaded with 0.8 mm OD, 5 mm long gold marker (1)
1325-3	18 gauge, 10 cm long Pre-waxed needles (3) preloaded with 0.8 mm OD, 7 mm long gold marker (1)
1325-4	18 gauge, 10 cm long Pre-waxed needles (3) preloaded with 0.8 mm OD, 10 mm long gold marker (1)
1325-5	18 gauge, 10 cm long Pre-waxed needles (3) preloaded with 0.8 mm OD, 1.5 mm long gold marker (1)
	Sterile Fiducial Marker Kit (1 per pack), 18 gauge, 20 cm Long Needle (0.8 mm OD)
1004-1	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 0.8 OD, 3 mm long gold marker (1)
1004-2	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 0.8 OD, 5 mm long gold marker (1)
1004-3	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 0.8 OD, 7 mm long gold marker (1)
1004-4	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 0.8 OD, 10 mm long gold marker (1)
	Sterile Fiducial Marker Kit (2 per pack), 18 gauge, 20 cm Long Needle (0.8 mm OD)
1430-1	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 0.8 OD, 3 mm long gold marker (1)
1430-2	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 0.8 OD, 5 mm long gold marker (1)
1430-3	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 0.8 OD, 7 mm long gold marker (1)
1430-4	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 0.8 OD, 10 mm long gold marker (1)
	Sterile Fiducial Marker Kit (3 per pack), 18 gauge, 20 cm Long Needle (0.8 mm OD)
1017-1	18 gauge, 20 cm long Pre-waxed needles (3) preloaded with 0.8 mm OD, 3 mm long gold marker (1)
1017-2	18 gauge, 20 cm long Pre-waxed needles (3) preloaded with 0.8 mm OD, 5 mm long gold marker (1)
1017-3	18 gauge, 20 cm long Pre-waxed needles (3) preloaded with 0.8 mm OD, 7 mm long gold marker (1)
1017-4	18 gauge, 20 cm long Pre-waxed needles (3) preloaded with 0.8 mm OD, 10 mm long gold markers (1)
	Sterile Fiducial Marker Kit (4 per pack), 18 gauge, 20 cm Long Needle (0.8 mm OD)
1168-1	18 gauge, 20 cm long Pre-waxed needle (4) preloaded with 0.8 OD, 3 mm long gold marker (1)
1168-2	18 gauge, 20 cm long Pre-waxed needle (4) preloaded with 0.8 OD, 5 mm long gold marker (1)
1168-3	18 gauge, 20 cm long Pre-waxed needle (4) preloaded with 0.8 OD, 7 mm long gold marker (1)
1168-4	18 gauge, 20 cm long Pre-waxed needle (4) preloaded with 0.8 OD, 10 mm long gold marker (1)

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
 phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



Part Number	Description
Sterile Fiducial Marker Kit (1 per pack), 18 gauge, 10 cm Long Needle (1.0 mm OD)	
1318-1	18 gauge, 10 cm long Pre-waxed needle (1) preloaded with 1.0 mm OD, 3 mm long gold marker (1)
1318-2	18 gauge, 10 cm long Pre-waxed needle (1) preloaded with 1.0 mm OD, 5 mm long gold marker (1)
1318-3	18 gauge, 10 cm long Pre-waxed needle (1) preloaded with 1.0 mm OD, 7 mm long gold marker (1)
1318-4	18 gauge, 10 cm long Pre-waxed needle (1) preloaded with 1.0 mm OD, 10 mm long gold marker (1)
Sterile Fiducial Marker Kit (2 per pack), 18 gauge, 10 cm Long Needle (1.0 mm OD)	
1432-1	18 gauge, 10 cm long Pre-waxed needle (2) preloaded with 1.0 mm OD, 3 mm long gold marker (1)
1432-2	18 gauge, 10 cm long Pre-waxed needle (2) preloaded with 1.0 mm OD, 5 mm long gold marker (1)
1432-3	18 gauge, 10 cm long Pre-waxed needle (2) preloaded with 1.0 mm OD, 7 mm long gold marker (1)
1432-4	18 gauge, 10 cm long Pre-waxed needle (2) preloaded with 1.0 mm OD, 10 mm long gold marker (1)
Sterile Fiducial Marker Kit (3 per pack), 18 gauge, 10 cm Long Needle (1.0 mm OD)	
1320-1	18 gauge, 10 cm long Pre-waxed needles (3) preloaded with 1.0 mm OD, 3 mm long gold marker (1)
1320-2	18 gauge, 10 cm long Pre-waxed needles (3) preloaded with 1.0 mm OD, 5 mm long gold marker (1)
1320-3	18 gauge, 10 cm long Pre-waxed needles (3) preloaded with 1.0 mm OD, 7 mm long gold marker (1)
1320-4	18 gauge, 10 cm long Pre-waxed needles (3) preloaded with 1.0 mm OD, 10 mm long gold marker (1)
Sterile Fiducial Marker Kit (1 per pack), 18 gauge, 20 cm Long Needle (1.0 mm OD)	
1008-1	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 1.0 mm OD, 3 mm long gold marker (1)
1008-2	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 1.0 mm OD, 5 mm long gold marker (1)
1008-3	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 1.0 mm OD, 7 mm long gold marker (1)
1008-4	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 1.0 mm OD, 10 mm long gold marker (1)
Sterile Fiducial Marker Kit (2 per pack), 18 gauge, 20 cm Long Needle (1.0 mm OD)	
1434-1	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 1.0 mm OD, 3 mm long gold marker (1)
1434-2	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 1.0 mm OD, 5 mm long gold marker (1)
1434-3	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 1.0 mm OD, 7 mm long gold marker (1)
1434-4	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 1.0 mm OD, 10 mm long gold marker (1)
Sterile Fiducial Marker Kit (3 per pack), 18 gauge, 20 cm Long Needle (1.0 mm OD)	
1019-1	18 gauge, 20 cm long Pre-waxed needles (3) preloaded with 1.0 mm OD, 3 mm long gold marker (1)
1019-2	18 gauge, 20 cm long Pre-waxed needles (3) preloaded with 1.0 mm OD, 5 mm long gold marker (1)
1019-3	18 gauge, 20 cm long Pre-waxed needles (3) preloaded with 1.0 mm OD, 7 mm long gold marker (1)
1019-4	18 gauge, 20 cm long Pre-waxed needles (3) preloaded with 1.0 mm OD, 10 mm long gold marker (1)

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
 phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



Part Number	Description
	Sterile Fiducial Marker Kit (4 per pack), 18 gauge, 20 cm Long Needle (1.0 mm OD)
1170-1	18 gauge, 20 cm long Pre-waxed needles (4) preloaded with 1.0 mm OD, 3 mm long gold marker (1)
1170-2	18 gauge, 20 cm long Pre-waxed needles (4) preloaded with 1.0 mm OD, 5 mm long gold marker (1)
1170-3	18 gauge, 20 cm long Pre-waxed needles (4) preloaded with 1.0 mm OD, 7 mm long gold marker (1)
1170-4	18 gauge, 20 cm long Pre-waxed needles (4) preloaded with 1.0 mm OD, 10 mm long gold marker (1)
	Sterile Fiducial Marker Kit (1 per pack), 17 gauge, 20 cm Long Needle (1.2 mm OD)
1088-1	17 gauge, 20 cm long Pre-waxed needles (1) preloaded with 1.2 OD, 3 mm long gold marker (1)
1088-2	17 gauge, 20 cm long Pre-waxed needles (1) preloaded with 1.2 OD, 5 mm long gold marker (1)
1088-3	17 gauge, 20 cm long Pre-waxed needles (1) preloaded with 1.2 OD, 7 mm long gold marker (1)
1088-4	17 gauge, 20 cm long Pre-waxed needles (1) preloaded with 1.2 OD, 10 mm long gold marker (1)
	Sterile Fiducial Marker Kit (2 per pack), 17 gauge, 20 cm Long Needle (1.2 mm OD)
1436-1	17 gauge, 20 cm long Pre-waxed needles (2) preloaded with 1.2 OD, 3 mm long gold marker (1)
1436-2	17 gauge, 20 cm long Pre-waxed needles (2) preloaded with 1.2 OD, 5 mm long gold marker (1)
1436-3	17 gauge, 20 cm long Pre-waxed needles (2) preloaded with 1.2 OD, 7 mm long gold marker (1)
1436-4	17 gauge, 20 cm long Pre-waxed needles (2) preloaded with 1.2 OD, 10 mm long gold marker (1)
	Sterile Fiducial Marker Kit (3 per pack), 17 gauge, 20 cm Long Needle (1.2 mm OD)
1090-1	17 gauge, 20 cm long Pre-waxed needles (3) preloaded with 1.2 OD, 3 mm long gold marker (1)
1090-2	17 gauge, 20 cm long Pre-waxed needles (3) preloaded with 1.2 OD, 5 mm long gold marker (1)
1090-3	17 gauge, 20 cm long Pre-waxed needles (3) preloaded with 1.2 OD, 7 mm long gold marker (1)
1090-4	17 gauge, 20 cm long Pre-waxed needles (3) preloaded with 1.2 OD, 10 mm long gold marker (1)
	Sterile Fiducial Marker Kit (4 per pack), 17 gauge, 20 cm Long Needle (1.2 mm OD)
1172-1	17 gauge, 20 cm long Pre-waxed needles (4) preloaded with 1.2 OD, 3 mm long gold marker (1)
1172-2	17 gauge, 20 cm long Pre-waxed needles (4) preloaded with 1.2 OD, 5 mm long gold marker (1)
1172-3	17 gauge, 20 cm long Pre-waxed needles (4) preloaded with 1.2 OD, 7 mm long gold marker (1)
1172-4	17 gauge, 20 cm long Pre-waxed needles (4) preloaded with 1.2 OD, 10 mm long gold marker (1)
	Sterile Fiducial Marker Kit (1 per pack), 17 gauge, 30 cm Long Needle (1.2 mm OD)
1031-1	17 gauge, 30 cm long Pre-waxed needles (1) preloaded with 1.2 OD, 3 mm long gold marker (1)
1031-2	17 gauge, 30 cm long Pre-waxed needles (1) preloaded with 1.2 OD, 5 mm long gold marker (1)
1031-3	17 gauge, 30 cm long Pre-waxed needles (1) preloaded with 1.2 OD, 7 mm long gold marker (1)
1031-4	17 gauge, 30 cm long Pre-waxed needles (1) preloaded with 1.2 OD, 10 mm long gold marker (1)

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
 phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



Part Number	Description
	Sterile Fiducial Marker Kit (2 per pack), 17 gauge, 30 cm Long Needle (1.2 mm OD)
1438-1	17 gauge, 30 cm long Pre-waxed needles (2) preloaded with 1.2 OD, 3 mm long gold marker (1)
1438-2	17 gauge, 30 cm long Pre-waxed needles (2) preloaded with 1.2 OD, 5 mm long gold marker (1)
1438-3	17 gauge, 30 cm long Pre-waxed needles (2) preloaded with 1.2 OD, 7 mm long gold marker (1)
1438-4	17 gauge, 30 cm long Pre-waxed needles (2) preloaded with 1.2 OD, 10 mm long gold marker (1)
	Sterile Fiducial Marker Kit (3 per pack), 17 gauge, 30 cm Long Needle (1.2 mm OD)
1033-1	17 gauge, 30 cm long Pre-waxed needles (3) preloaded with 1.2 OD, 3 mm long gold marker (1)
1033-2	17 gauge, 30 cm long Pre-waxed needles (3) preloaded with 1.2 OD, 5 mm long gold marker (1)
1033-3	17 gauge, 30 cm long Pre-waxed needles (3) preloaded with 1.2 OD, 7 mm long gold marker (1)
1033-4	17 gauge, 30 cm long Pre-waxed needles (3) preloaded with 1.2 OD, 10 mm long gold marker (1)
	Sterile Radiopaque Marker Strand (1 per pack), 18 gauge, 10 cm Long Needle (0.8 mm OD)
1293-1-X	18 gauge, 10 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 0.8 mm OD, 3 mm long gold marker
1293-2-X	18 gauge, 10 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 0.8 mm OD, 5 mm long gold marker
1293-3-X	18 gauge, 10 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 0.8 mm OD, 7 mm long gold marker
1293-4-X	18 gauge, 10 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 0.8 mm OD, 10 mm long gold marker
	Sterile Radiopaque Marker Strand (2 per pack), 18 gauge, 10 cm Long Needle (0.8 mm OD)
1440-1-X	18 gauge, 10 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 0.8 mm OD, 3 mm long gold marker
1440-2-X	18 gauge, 10 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 0.8 mm OD, 5 mm long gold marker
1440-3-X	18 gauge, 10 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 0.8 mm OD, 7 mm long gold marker
1440-4-X	18 gauge, 10 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 0.8 mm OD, 10 mm long gold marker
	Sterile Radiopaque Marker Strand (3 per pack), 18 gauge, 10 cm Long Needle (0.8 mm OD)
1295-1-X	18 gauge, 10 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 0.8 mm OD, 3 mm long gold marker
1295-2-X	18 gauge, 10 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 0.8 mm OD, 5 mm long gold marker
1295-3-X	18 gauge, 10 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 0.8 mm OD, 7 mm long gold marker
1295-4-X	18 gauge, 10 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 0.8 mm OD, 10 mm long gold marker

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
 phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



Part Number	Description
	Sterile Radiopaque Marker Strand (1 per pack), 18 gauge, 20 cm Long Needle (0.8 mm OD)
1013-1-X	18 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 0.8 mm OD, 3 mm long gold marker
1013-2-X	18 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 0.8 mm OD, 5 mm long gold marker
1013-3-X	18 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 0.8 mm OD, 7 mm long gold marker
1013-4-X	18 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 0.8 mm OD, 10 mm long gold marker
	Sterile Radiopaque Marker Strand (2 per pack), 18 gauge, 20 cm Long Needle (0.8 mm OD)
1442-1-X	18 gauge, 20 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 0.8 mm OD, 3 mm long gold marker
1442-2-X	18 gauge, 20 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 0.8 mm OD, 5 mm long gold marker
1442-3-X	18 gauge, 20 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 0.8 mm OD, 7 mm long gold marker
1442-4-X	18 gauge, 20 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 0.8 mm OD, 10 mm long gold marker
	Sterile Radiopaque Marker Strand (3 per pack), 18 gauge, 20 cm Long Needle (0.8 mm OD)
1015-1-X	18 gauge, 20 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 0.8 mm OD, 3 mm long gold marker
1015-2-X	18 gauge, 20 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 0.8 mm OD, 5 mm long gold marker
1015-3-X	17 gauge, 20 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 0.8 mm OD, 7 mm long gold marker
1015-4-X	17 gauge, 20 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 0.8 mm OD, 10 mm long gold marker
	Sterile Radiopaque Marker Strand (1 per pack), 17 gauge, 30 cm Long Needle (0.8 mm OD)
1130-1-X	17 gauge, 30 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 0.8 mm OD, 3 mm long gold marker
1130-2-X	17 gauge, 30 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 0.8 mm OD, 5 mm long gold marker
1130-3-X	17 gauge, 30 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 0.8 mm OD, 7 mm long gold marker
1130-4-X	17 gauge, 30 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 0.8 mm OD, 10 mm long gold marker
	Sterile Radiopaque Marker Strand (2 per pack), 17 gauge, 30 cm Long Needle (0.8 mm OD)
1444-1-X	17 gauge, 30 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 0.8 mm OD, 3 mm long gold marker
1444-2-X	17 gauge, 30 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 0.8 mm OD, 5 mm long gold marker

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
 phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



Part Number	Description
1444-3-X	17 gauge, 30 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 0.8 mm OD, 7 mm long gold marker
1444-4-X	17 gauge, 30 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 0.8 mm OD, 10 mm long gold marker
	Sterile Radiopaque Marker Strand (3 per pack), 17 gauge, 30 cm Long Needle (0.8 mm OD)
1132-1-X	17 gauge, 30 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 0.8 mm OD, 3 mm long gold marker
1132-2-X	17 gauge, 30 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 0.8 mm OD, 5 mm long gold marker
1132-3-X	17 gauge, 30 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 0.8 mm OD, 7 mm long gold marker
1132-4-X	17 gauge, 30 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 0.8 mm OD, 10 mm long gold marker
	Sterile Radiopaque Marker Strand (1 per pack), 17 gauge, 10 cm Long Needle (1.0 mm OD)
1357-1-X	17 gauge, 10 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.0 mm OD, 3 mm long gold marker
1357-2-X	17 gauge, 10 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.0 mm OD, 5 mm long gold marker
1357-3-X	17 gauge, 10 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.0 mm OD, 7 mm long gold marker
1357-4-X	17 gauge, 10 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.0 mm OD, 10 mm long gold marker
1357-5-X	17 gauge, 10 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.0 mm OD, 1.5 mm long gold marker
	Sterile Radiopaque Marker Strand (2 per pack), 17 gauge, 10 cm Long Needle (1.0 mm OD)
1446-1-X	17 gauge, 10 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.0 mm OD, 3 mm long gold marker
1446-2-X	17 gauge, 10 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.0 mm OD, 5 mm long gold marker
1446-3-X	17 gauge, 10 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.0 mm OD, 7 mm long gold marker
1446-4-X	17 gauge, 10 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.0 mm OD, 10 mm long gold marker
1446-5-X	17 gauge, 10 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.0 mm OD, 1.5 mm long gold marker
	Sterile Radiopaque Marker Strand (3 per pack), 17 gauge, 10 cm Long Needle (1.0 mm OD)
1359-1-X	17 gauge, 10 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.0 mm OD, 3 mm long gold marker
1359-2-X	17 gauge, 10 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.0 mm OD, 5 mm long gold marker

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
 phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



Part Number	Description
1359-3-X	17 gauge, 10 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.0 mm OD, 7 mm long gold marker
1359-4-X	17 gauge, 10 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.0 mm OD, 10 mm long gold marker
1359-5-X	17 gauge, 10 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.0 mm OD, 1.5 mm long gold marker
	Sterile Radiopaque Marker Strand (1 per pack), 17 gauge, 20 cm Long Needle (1.0 mm OD)
1268-1-X	17 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.0 mm OD, 3 mm long gold marker
1268-2-X	17 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.0 mm OD, 5 mm long gold marker
1268-3-X	17 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.0 mm OD, 7 mm long gold marker
1268-4-X	17 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.0 mm OD, 10 mm long gold marker
1268-5-X	17 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.0 mm OD, 1.5 mm long gold marker
	Sterile Radiopaque Marker Strand (2 per pack), 17 gauge, 20 cm Long Needle (1.0 mm OD)
1448-1-X	17 gauge, 20 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.0 mm OD, 3 mm long gold marker
1448-2-X	17 gauge, 20 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.0 mm OD, 5 mm long gold marker
1448-3-X	17 gauge, 20 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.0 mm OD, 7 mm long gold marker
1448-4-X	17 gauge, 20 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.0 mm OD, 10 mm long gold marker
1448-5-X	17 gauge, 20 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.0 mm OD, 1.5 mm long gold marker
	Sterile Radiopaque Marker Strand (3 per pack), 17 gauge, 20 cm Long Needle (1.0 mm OD)
1270-1-X	17 gauge, 20 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.0 mm OD, 3 mm long gold marker
1270-2-X	17 gauge, 20 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.0 mm OD, 5 mm long gold marker
1270-3-X	17 gauge, 20 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.0 mm OD, 7 mm long gold marker
1270-4-X	17 gauge, 20 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.0 mm OD, 10 mm long gold marker
1270-5-X	17 gauge, 20 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.0 mm OD, 1.5 mm long gold marker
	Sterile Radiopaque Marker Strand (1 per pack), 16 gauge, 20 cm Long Needle (1.2 mm OD)
1329-1-X	16 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.2 mm OD, 3 mm long gold marker

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
 phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



Part Number	Description
1329-2-X	16 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.2 mm OD, 5 mm long gold marker
1329-3-X	16 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.2 mm OD, 7 mm long gold marker
1329-4-X	16 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.2 mm OD, 10 mm long gold marker
1329-5-X	16 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 1.2 mm OD, 1.5 mm long gold marker
	Sterile Radiopaque Marker Strand (2 per pack), 16 gauge, 20cm Long Needle (1.2 mm OD)
1450-1-X	16 gauge, 20 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.2 mm OD, 3 mm long gold marker
1450-2-X	16 gauge, 20 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.2 mm OD, 5 mm long gold marker
1450-3-X	16 gauge, 20 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.2 mm OD, 7 mm long gold marker
1450-4-X	16 gauge, 20 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.2 mm OD, 10 mm long gold marker
1450-5-X	16 gauge, 20 cm long Pre-waxed needles (2) preloaded with absorbable strand enclosing 1.2 mm OD, 1.5 mm long gold marker
	Sterile Radiopaque Marker Strand (3 per pack), 16 gauge, 20cm Long Needle (1.2 mm OD)
1331-1-X	16 gauge, 20 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.2mm OD, 3 mm long gold marker
1331-2-X	16 gauge, 20 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.2 mm OD, 5 mm long gold marker
1331-3-X	16 gauge, 20 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.2 mm OD, 7 mm long gold marker
1331-4-X	16 gauge, 20 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.2 mm OD, 10 mm long gold marker
1331-5-X	16 gauge, 20 cm long Pre-waxed needles (3) preloaded with absorbable strand enclosing 1.2 mm OD, 1.5 mm long gold marker
	Sterile Needle w/Best Marker Strand w/Iodine Dummy Seed (1 per pack), 18 gauge, 20 cm Long Needle (0.8 mm OD)
10113-17-X	18 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 0.8 mm OD, 5 mm long iodine dummy seed
	Sterile Needle w/Best Marker Strand w/Iodine Dummy Seed (3 per pack), 18 gauge, 20 cm Long Needle (0.8 mm OD)
10113-19-X	18 gauge, 20 cm long Pre-waxed needles (1) preloaded with absorbable strand enclosing 0.8 mm OD, 5 mm long iodine dummy seed
	Sterile Fiducial Marker Kit With Spacers (1 per pack), 18 gauge, 20 cm Long Needle (0.8 mm OD)
1306-1-X	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 0.8 OD, 3 mm long gold marker(s) and Spacer(s)

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
 phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



Part Number	Description
1306-2-X	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 0.8 OD, 5 mm long gold marker(s) and Spacer(s)
1306-3-X	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 0.8 OD, 7 mm long gold marker(s) and Spacer(s)
1306-4-X	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 0.8 OD, 10 mm long gold marker(s) and Spacer(s)
Sterile Fiducial Marker Kit With Spacers (2 per pack), 18 gauge, 20 cm Long Needle (0.8 mm OD)	
1452-1-X	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 0.8 OD, 3 mm long gold marker(s) and Spacer(s)
1452-2-X	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 0.8 OD, 5 mm long gold marker(s) and Spacer(s)
1452-3-X	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 0.8 OD, 7 mm long gold marker(s) and Spacer(s)
1452-4-X	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 0.8 OD, 10 mm long gold marker(s) and Spacer(s)
Sterile Fiducial Marker Kit With Spacers (3 per pack), 18 gauge, 20 cm Long Needle (0.8 mm OD)	
1385-1-X	18 gauge, 20 cm long Pre-waxed needle (3) preloaded with 0.8 OD, 3 mm long gold marker(s) and Spacer(s)
1385-2-X	18 gauge, 20 cm long Pre-waxed needle (3) preloaded with 0.8 OD, 5 mm long gold marker(s) and Spacer(s)
1385-3-X	18 gauge, 20 cm long Pre-waxed needle (3) preloaded with 0.8 OD, 7 mm long gold marker(s) and Spacer(s)
1385-4-X	18 gauge, 20 cm long Pre-waxed needle (3) preloaded with 0.8 OD, 10 mm long gold marker(s) and Spacer(s)
Sterile Fiducial Marker Kit With Spacers (1 per pack), 18 gauge, 20 cm Long Needle (1.0 mm OD)	
1309-1-X	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 1.0 OD, 3 mm long gold marker(s) and Spacer(s)
1309-2-X	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 1.0 OD, 5 mm long gold marker(s) and Spacer(s)
1309-3-X	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 1.0 OD, 7 mm long gold marker(s) and Spacer(s)
1309-4-X	18 gauge, 20 cm long Pre-waxed needle (1) preloaded with 1.0 OD, 10 mm long gold marker(s) and Spacer(s)
Sterile Fiducial Marker Kit With Spacers (2 per pack), 18 gauge, 20 cm Long Needle (1.0 mm OD)	
1454-1-X	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 1.0 OD, 3 mm long gold marker(s) and Spacer(s)
1454-2-X	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 1.0 OD, 5 mm long gold marker(s) and Spacer(s)
1454-3-X	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 1.0 OD, 7 mm long gold marker(s) and Spacer(s)

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
 phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



Part Number	Description
1454-4-X	18 gauge, 20 cm long Pre-waxed needle (2) preloaded with 1.0 OD, 10 mm long gold marker(s) and Spacer(s)
	Sterile Fiducial Marker Kit With Spacers (3 per pack), 18 gauge, 20 cm Long Needle (1.0 mm OD)
1387-1-X	18 gauge, 20 cm long Pre-waxed needle (3) preloaded with 1.0 OD, 3 mm long gold marker(s) and Spacer(s)
1387-2-X	18 gauge, 20 cm long Pre-waxed needle (3) preloaded with 1.0 OD, 5 mm long gold marker(s) and Spacer(s)
1387-3-X	18 gauge, 20 cm long Pre-waxed needle (3) preloaded with 1.0 OD, 7 mm long gold marker(s) and Spacer(s)
1387-4-X	18 gauge, 20 cm long Pre-waxed needle (3) preloaded with 1.0 OD, 10 mm long gold marker(s) and Spacer(s)
	Sterile Fiducial Marker Kit With Spacers (1 per pack), 17 gauge, 20 cm Long Needle (1.2 mm OD)
1312-1-X	17 gauge, 20 cm long Pre-waxed needle (1) preloaded with 1.2 OD, 3 mm long gold marker(s) and Spacer(s)
1312-2-X	17 gauge, 20 cm long Pre-waxed needle (1) preloaded with 1.2 OD, 5 mm long gold marker(s) and Spacer(s)
1312-3-X	17 gauge, 20 cm long Pre-waxed needle (1) preloaded with 1.2 OD, 7 mm long gold marker(s) and Spacer(s)
1312-4-X	17 gauge, 20 cm long Pre-waxed needle (1) preloaded with 1.2 OD, 10 mm long gold marker(s) and Spacer(s)
	Sterile Fiducial Marker Kit With Spacers (2 per pack), 17 gauge, 20 cm Long Needle (1.2 mm OD)
1456-1-X	17 gauge, 20 cm long Pre-waxed needle (2) preloaded with 1.2 OD, 3 mm long gold marker(s) and Spacer(s)
1456-2-X	17 gauge, 20 cm long Pre-waxed needle (2) preloaded with 1.2 OD, 5 mm long gold marker(s) and Spacer(s)
1456-3-X	17 gauge, 20 cm long Pre-waxed needle (2) preloaded with 1.2 OD, 7 mm long gold marker(s) and Spacer(s)
1456-4-X	17 gauge, 20 cm long Pre-waxed needle (2) preloaded with 1.2 OD, 10 mm long gold marker(s) and Spacer(s)
	Sterile Fiducial Marker Kit With Spacers (3 per pack), 17 gauge, 20 cm Long Needle (1.2 mm OD)
1389-1-X	17 gauge, 20 cm long Pre-waxed needle (3) preloaded with 1.2 OD, 3 mm long gold marker(s) and Spacer(s)
1389-2-X	17 gauge, 20 cm long Pre-waxed needle (3) preloaded with 1.2 OD, 5 mm long gold marker(s) and Spacer(s)
1389-3-X	17 gauge, 20 cm long Pre-waxed needle (3) preloaded with 1.2 OD, 7 mm long gold marker(s) and Spacer(s)
1389-4-X	17 gauge, 20 cm long Pre-waxed needle (3) preloaded with 1.2 OD, 10 mm long gold marker(s) and Spacer(s)

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
 phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



View

Best Medical International's
Complete Catalog:

www.bestmedical.com/products/catalog.pdf