

Iodine-125 and Palladium-103

Patented Double Wall Best® Iodine-125 Source

Model # 2301



Best® Iodine-125 and Best® Palladium-103 Seeds:

- Double Wall Encapsulation
- Traceable to NIST₁₉₉₉ Calibration Revised for 2000
- Sterile Absorbable Strands (Customer's Prescription)
- Sterile, Pre-Loaded Needles (Customer's Choice of Needles)
- Loaded in Mick Cartridges, Sterile or Non-Sterile
- Loose Seeds

5 mm Seeds — 5 mm Spacers

- Seeds and spacers are the same length (5mm). Double wall encapsulation and uniform wall thickness with superior dose distribution. Consistent quality, and seeds fit Mick cartridges.
- Excellent visualization with Ultrasound, X-ray, Fluoro, CT and MRI. **Real time dosimetry possible with Best® seeds.**
- Complies with AAPM's dosimetric prerequisites

Patented Double Wall Best® Palladium-103 Source



Model # 2335

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com


Best[®]
healthcare for everyone

Best® Iodine-125 and Best® Palladium-103 Sources

- Patented Double Wall Titanium Encapsulation with uniform wall thickness
- Tungsten Marker—excellent images of source in ultrasound, X-ray, CT, MRI, and fluoro—far superior to other commercially available sources
- Source Calibration traceable to NIST₁₉₉₉ Calibration revised for 2000
- Sources and spacers are the same length (5.0 mm nominal) and diameter (0.8 mm nominal), fitting Mick cartridge snugly
- Custom made high quality spacers and Mick cartridges provided at no cost
- Sources can be provided as loose seeds; loaded in needles with spacers, sterile; loaded in strands, sterile; preloaded in Mick cartridges

Patented Double Wall Best® Iodine-125 Source (Model 2301)

Source Trade Name: **Best® Iodine Source**
Marker: **Tungsten**
Model #: **2301**
Dose Rate Constant: **(NIST 1999 with 2000 modification)**
External Length: **5 mm (Nominal) 1.02 cGy/hr/U**
External Diameter: **0.8 mm (Nominal)**

Patented Double Wall Best® Palladium-103 Source (Model 2335)

Source Trade Name: **Best® Palladium Source**
Marker: **Tungsten**
Model #: **2335**
Dose Rate Constant: **(NIST 1999 with 2000 modification)**
External Length: **5 mm (Nominal) 0.70 cGy/hr/U**
External Diameter: **0.8 mm (Nominal)**

To Place an Order, Call (800) 336-4970 or Fax (703) 451-4852

References:

- 1) A.S. Meigooni, D.M. Gearheart, and K. Sowards, Med. Phys. 27, 2168-2173 (2000)
- 2) S.W. Peterson and B. Thomadsen, Brachytherapy 1, 110-119 (2002)
- 3) R. Nath and N. Yue, Brachytherapy, 1, 102-109 (2002)

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



Sterile Seed Products

I-125

Part #	Description
10165-1	I-125 Strands in a Pack
10165-3	Needle w/I-125 Sources in Strands in a Tray
10165-4	Needle w/I-125 Loose Sources in a Tray
10165-17	Needle w/I-125 Loose Sources and Strands in a Tray
10113-11	Cartridge with I-125 Sources in a Pouch



Pd-103

Part #	Description
15165-1	Pd-103 Strands in a Pack
15165-4	Needle w/Pd-103 Loose Sources in a Tray
15165-9	Needle w/Pd-103 Sources in Strands in a Tray
15165-11	Needle w/Pd-103 Loose Sources and Strands in a Tray
15113-4	Cartridge with Pd-103 Sources in a Pouch



Available from Best® within 24 hours:
Best® I-125 and Pd-103 seeds in sterile packing



- Best® Custom Strand (Per Your Prescription)
- Best® Standard Strand (10 Seeds/Strand Spaced 1 cm Apart)
- Best® Real Time Strand/Needle (Use Our Loader in the OR)
- Best® Seeds and Spacers in Pre-loaded Needles
- Cartridges
- Loose Seeds

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



Best ADVANTAGE

- We are the largest supplier of Iridium 192 and Gold 198 seeds in the world.
- We are the first and only provider of Iridium 192 in nylon ribbons in a wide range of activities from 0.10 to 10 mg Raeq per seed or 0.3 to 30 mg Raeq per centimeter.
- We are the first and only company to provide custom service for Gold 198 permanent implants as well as a wide range of activities from 0.1 to 30 mCi per seed.
- We are the only company to provide a number of custom-made accessories for temporary and permanent implants in Brachytherapy.
- We are the first and only company to call our customers the morning after shipping to inform them of mode of transportation and expected arrival time.
- We can supply Ir-192, Au-198, Pd-103, I-125 and Brachytherapy accessories on a single Purchase Order.
- We ship most materials within 24 hours of receipt of the order.

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
phone 703 451 2378 800 336 4970 fax 703 451 5228 www.bestmedical.com



View

Best Medical International's
Complete Catalog:

www.bestmedical.com/products/catalog.pdf